



Amtryke Adaptive Tricycle Request Form

(Must be filled out completely by adult rider or parent/guardian)

Recipient's Name: _____ Age: _____ Date of Request: _____
Mailing Address: _____ Phone #: _____
City/State/Zip: _____ Email: _____
Diagnosis: _____

If Recipient is Under Age 18

Parent/Guardian Name: _____
If different from above
Mailing Address: _____ Phone #: _____
City/State/Zip: _____ Email: _____

Secondary Contact Name: _____ Phone #: _____

Treating Therapist's Name: _____

Phone #: _____ Email: _____

How did you hear about the Amtryke Adaptive Tricycle? (Check all that apply)

Therapist Website AMBUCS Member Other: _____

Tell us about the recipient*:

Signature: _____ Date: _____

Before an order can be placed, we must receive a fully completed Request & Waiver Form, Assessment Form and Model Selection Form.

Please e-mail completed forms to caprocktrykes@gmail.com or mail copy to P.O. Box 93871, Lubbock, Tx 79423.

You will receive either a confirmation e-mail from caprocktrykes@gmail.com or a phone call from a member to inform you that your application has been received.

Amtryke Adaptive Tricycle Waiver Form

(Must be filled out completely by adult rider or parent/guardian)

AMBUCS members nationwide are dedicated to creating opportunities for mobility and independence by providing Amtryke adaptive tricycles, offering educational scholarships to therapy students and performing various forms of community service.

Purpose: The Amtryke adaptive tricycle creates a feeling of freedom, builds self-esteem, strengthens muscles and improves motor coordination and range of motion—all while making exercise fun.

Steering: Initially, the rider may have difficulty turning or changing directions. Encourage the rider to go straight ahead, back up and slowly turn around. On many models there are three steering options for the Amtryke. On the front column of the tricycle you will find two holes for the steering pin. The top hole is straight steering, the bottom hole allows a 20-degree turning radius. Leaving the pin out gives the rider free steering.

Safety Cautions

- Fast speeds and sharp turns can cause the Amtryke adaptive tricycle to tip or turn over.
- Always wear a helmet when riding an Amtryke. Use of other protective gear is highly recommended.
- Adult supervision required if used by younger or developmentally delayed riders.
- Use caution near vehicles, swimming pools and other bodies of water, hills, alleys and sloped driveways.
- Always wear shoes.
- Never allow more than one rider.
- Use of the steering pin is recommended to prevent over-steering or possible tip-overs.

The information contained in this service is not intended nor implied by National AMBUCS™, Inc. to be professional medical advice by National AMBUCS, Inc. Always seek the advice of your physician, therapist or other qualified healthcare provider prior to starting any treatment or with any question you may have regarding a medical condition. Nothing contained in this document is intended by National AMBUCS, Inc. to be for medical diagnosis or treatment by National AMBUCS, Inc. or on behalf of National AMBUCS, Inc.

In no event shall National AMBUCS, Inc. be liable for any direct, indirect, incidental, consequential, special, exemplary, punitive, or any other monetary or other damages, personal injury or property damages, fees, fines, costs, attorney fees, or liabilities of any kind arising out of or relating in any way to this service or use of the Amtryke® adaptive tricycle, and/or content or information provided herein.

I agree that the rider listed on the front of this document may be photographed. I also agree that the photos and recipient's name may be used in promotional efforts for National AMBUCS, Amtryke or the local AMBUCS chapter. I further grant AMBUCS the ability to use the photos and name for advertising/publicity purposes without additional compensation, except where prohibited by law. If anybody in my party does not want to be photographed under these same terms, I will let the photographer know as soon as possible.

By signing below, I acknowledge that I have read and understood this liability waiver.

Recipient's Name: _____

Adult Recipient Signature: _____

If Recipient is Under Age 18

Legal Guardian Name: _____

Legal Guardian Signature: _____ Date: _____

Before an order can be placed, we must receive a fully completed Request & Waiver Form, Assessment Form and Model Selection Form.

Please e-mail completed forms to caprocktrykes@gmail.com or mail copy to P.O. Box 93871, Lubbock, Tx 79423.

You will receive either a confirmation e-mail from caprocktrykes@gmail.com or a phone call from a member to inform you that your application has been received.


Amtryke Assessment Form

(Must be filled out completely by therapist)

Recipient's Name: _____
 Age: _____ Weight (lbs.):* _____ Height (inches):* _____
 Diagnosis(es):* _____

*This private information is only used to help appropriately fit the rider.

RIDER'S MEASUREMENTS				
Arm Measurements (inches)				Total Length
Left	A to B:		B to C:	
Right	A to B:		B to C:	
Trunk	A to D:			
Leg Measurements (inches)				Total Length
Left	D to E:		E to F:	
Right	D to E:		E to F:	



A	Center of Shoulder
B	Center of Elbow
C	Center of Digit Crease
D	Center of Hip
E	Center of Knee
F	Bottom of Foot

Arm Length & Leg Length Measurements are critical to correct Amtryke Selection

Sizing Chart is available online:
www.ambucs.org/riders/wish-list/sizing-chart/

Notes on Provided Measurements (if any):

Helmet Sizing	
Size	Measurement (head circumference)
XXS	18.5" to 19.5"
XS	20.5" to 22"
S/M	22" to 23.6"
L/XL	23.6" to 25.75"

Therapist Name: _____ Is this the treating therapist? Yes No

Credentials: _____

Phone: _____ Email: _____

Facility Name: _____

Street Address: _____ City: _____ State: _____ Zip: _____

Is this facility an Amtryke Evaluation Site? Yes No Not Sure

Therapist comments concerning recipient or goals:

This request/assessment is directed to:

- Local AMBUCS Chapter Name: _____
 National Wish List (AMBUCS Resource Center)

By signing below, you are signifying that in your professional opinion this rider would benefit from an Amtryke. You assume no liability.

Therapist Signature: _____ **Date:** _____

Before an order can be placed, we must receive a fully completed Request & Waiver Form, Assessment Form and Model Selection Form.

Please e-mail completed forms to caprocktrykes@gmail.com or mail copy to P.O. Box 93871, Lubbock, Tx 79423.

You will receive either a confirmation e-mail from caprocktrykes@gmail.com or a phone call from a member to inform you that your application has been received.

In order to accommodate the widest variety of people, Amtryke offers many tryke models and each can be customized in a variety of ways. Following the steps below will help you choose to the perfect tryke for your client from what might seem like a dizzying array of options.

Remember you can always refer to our website, www.amtrykestore.org, or the Amtryke catalogue for more information and product images.

Step 1: Fill out the first page of the Amtryke Assessment Form.

Step 2: Choose the way the tryke will be propelled: **Hand & Foot, Foot, or Hand**. Your choice should be based on the rider's ability and therapy goals.

Hand & Foot trykes improve coordination, strength and range of motion. Using all four extremities helps with weakness in any area, even general weakness, and can positively affect tone.

Foot trykes were developed in response to requests from therapists for a traditional tricycle for riders with special needs.

Hand trykes are designed for persons whose lower limbs lack function or those who need special therapy for the upper extremity.

Step 3: Take rider's measurements from the front of this Assessment Form and compare them to our Sizing Chart. (This will narrow the choices considerably)

Step 4: Choose a drive. (If it doesn't mention a choice, then ignore this step.)

Amtrykes come with two drive possibilities: **fixed** drive or **geared** drive. Tryke models have been carefully designed so the drive matches the therapeutic goals of the equipment; therefore all drives are not available on all trykes.

A **fixed** drive, commonly known as a 'fixie,' works on a mechanical level to help individuals make a full pedal rotation. The foot crank is constantly in motion for full therapeutic effect. Coasting is not possible with a positive drive; when limb motion stops, the bike does as well.

Riders of **geared** drive trykes must be able to make a complete pedal rotation on their own. This tryke is suitable for riders who need help with balance and a stable sitting position and who have the cognitive ability to steer successfully and apply the coaster brake or hand brake. The key feature of geared trykes is the ability to coast.

Step 5: Chose any adaptations and/or accessories needed by the rider. *Each tryke model can only be customized in the ways noted in its own model section in the Tryke Selection Forms or with the generic accessories listed below. If a customization option is not listed, it is because of design or other practical constraints in offering it on a particular model.*

GENERIC ACCESSORIES (not model specific)			
Fun Items	<input type="checkbox"/> License Plate	<input type="checkbox"/> Water Bottle w Cage	<input type="checkbox"/> Mirror
Leg & Foot Items	<input type="checkbox"/> Foot Cups (pair): <input type="checkbox"/> Small <input type="checkbox"/> Medium	<input type="checkbox"/> Knee Adductor Strap: <input type="checkbox"/> Small <input type="checkbox"/> Large <input type="checkbox"/> Medium	<input type="checkbox"/> Pedal Block (1 = ¾") _____ qty
Hand Items	<input type="checkbox"/> Variable Range of Motion Kit <i>(only for Hand & Foot Cycles)</i>	<input type="checkbox"/> Wrist Wraps <i>(Includes right & left)</i> <input type="checkbox"/> X-Small <input type="checkbox"/> Large <input type="checkbox"/> Small <input type="checkbox"/> Medium	<input type="checkbox"/> Wrist Brace Mitt: <input type="checkbox"/> Right <input type="checkbox"/> Left <input type="checkbox"/> XX-Small <input type="checkbox"/> Medium <input type="checkbox"/> X-Small <input type="checkbox"/> Small

Before an order can be placed, we must receive a fully completed Request & Waiver Form, Assessment Form and Model Selection Form.

Please e-mail completed forms to caprocktrykes@gmail.com or mail copy to P.O. Box 93871, Lubbock, Tx 79423.

You will receive either a confirmation e-mail from caprocktrykes@gmail.com or a phone call from a member to inform you that your application has been received.

Recipient's Name: _____

FOOT TRYKES

<input type="checkbox"/> 1410 FOOT TRYKE – (10" wheels, fixed drive) 50-FC-0100 Arm Length 12-20", Leg Length 15-21" NEW: Standard with Rear Steering (can convert to push bar) and safety flag
<input type="checkbox"/> Option #1 – Standard Seating System: Blue Bucket Seat
<input type="checkbox"/> Option #2 – Alternate Seating System: Snappy Seat System Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest <input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back**
<input type="checkbox"/> Option #3 – Alternate Seating System: OLD Blue Bucket Seat** <i>(Less expensive option for very small riders. Cannot be used with an H-Harness or Separator Cube.)</i>
1410 Accessories: <input type="checkbox"/> Separator Cube <input type="checkbox"/> H-Harness - 11.5" <input type="checkbox"/> Heavy-Duty Push Bar – <i>for use with notched seats</i>

<input type="checkbox"/> 1412 ProSeries FOOT TRYKE (12" wheels, fixed drive) – 50-FC-1412 Arm Length 16-24", Leg Length 23-28" NEW: Standard with Rear Steering (can convert to push bar), Medium Pommel Saddle Seat and safety flag
<input type="checkbox"/> Option #1 – Standard Seating System: Medium Pommel Saddle Seat/1600 Simple Seatback Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Large Pommel Saddle Seat <input type="checkbox"/> Saddle Seat <input type="checkbox"/> Skinny Saddle Seat (8.5x6")
<input type="checkbox"/> Option #2 – Alternate Seating System: 1400 ProSeries Seat Back System (push grip, backpad, med. pommel & 2 laterals) Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Large Pommel Saddle Seat <input type="checkbox"/> Saddle Seat <input type="checkbox"/> Skinny Saddle Seat (8.5x6") Trunk Support Options: <input type="checkbox"/> ProSeries Full Padded Back <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Lumbar Pad <input type="checkbox"/> Recumbent 10 Degree Seat Post
<input type="checkbox"/> Option #3 – Alternate Seating System: Gray Bucket Seat
<input type="checkbox"/> Option #4 – Alternate Seating System: Snappy Seat System <i>(too small for many 1412 riders; check seat measurements.)</i> Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest <input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back** <input type="checkbox"/> Separator Cube
1412 Accessories: <input type="checkbox"/> Heavy-Duty Push Bar <input type="checkbox"/> Pedal Toe Pulley <input type="checkbox"/> H-Harness - 11.5" <input type="checkbox"/> Solid Rear Tires <input type="checkbox"/> Front Hand Brake Kit (cannot be used with Rear Steering) <input type="checkbox"/> ½" Expanding Pedals <input type="checkbox"/> Hi-Rise Handlebars: <input type="checkbox"/> 5" (BMX) <input type="checkbox"/> 8" <input type="checkbox"/> ½" Exercise Pedals <input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5" <input type="checkbox"/> 9" Extender Tube

*cannot be used in conjunction with rear steering. **cannot be used in conjunction with an h-harness, head rest or laterals.

<input type="checkbox"/> NEW! High Functioning Set Up (Both Options Allow Coasting) Can be used with 1416, 1420, 1420XL (Rider must have the ability to steer, apply brake and make a complete pedal rotation on their own.) Comes Standard with Pommel Saddle Seat (1412 – M; 1416 & 1420 – L), 1600 Simple Seat Back, Hi-Rise Handlebars, Exercise pedals and safety flag. MUST choose Option 1 or 2.
Choose Tryke: <input type="checkbox"/> 1416 <input type="checkbox"/> 1420 <input type="checkbox"/> 1420XL
<input type="checkbox"/> Option #1: (switch to existing freewheel hub and add Disk Brake Kit - <i>necessary</i> for safety) Adds the ability to coast <input type="checkbox"/> Dual Hand Brake (Option 1 only)
<input type="checkbox"/> Option #2: (3-Speed Kit and coaster brake) Adds the ability to coast and switch gears for hilly terrain – <i>Replaces the 1600 Series</i>
Standard Seating System (for either option): Pommel Saddle Seat/1600 Simple Seat Back Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Tractor Seat with bracket <input type="checkbox"/> Regular Saddle (11x9") <input type="checkbox"/> Medium Saddle (10.5x10") <input type="checkbox"/> Large Saddle (13x12")
Adaptive Accessories: <input type="checkbox"/> ½" Expanding Pedals <input type="checkbox"/> Telescoping Loop Handlebars <input type="checkbox"/> 1400 Seat Back Set (w/ laterals)

Before an order can be placed, we must receive a fully completed Request & Waiver Form, Assessment Form and Model Selection Form.

Please e-mail completed forms to caprocktrykes@gmail.com or mail copy to P.O. Box 93871, Lubbock, Tx 79423.

You will receive either a confirmation e-mail from caprocktrykes@gmail.com or a phone call from a member to inform you that your application has been received.

Recipient's Name:

FOOT TRYKES

<input type="checkbox"/> 1416 ProSeries FOOT TRYKE (16" wheels, fixed or freewheel drive) – 50-FC-1416 Arm Length 16-24", Leg Length 26-32" NEW: Standard with Medium Pommel Saddle Seat and safety flag		
<input type="checkbox"/> Option #1 – Standard Seating System: 1400 ProSeries Seat Back System (push grip, backpad, med. pommel & 2 laterals) Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Tractor Seat with bracket <input type="checkbox"/> Large Pommel Saddle Seat <input type="checkbox"/> Saddle Seat Trunk Support Options: <input type="checkbox"/> ProSeries Full Padded Back <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Lumbar Pad <input type="checkbox"/> Recumbent 10 Degree Seat Post		
<input type="checkbox"/> Option #2 – Alternate Seating System: Medium Pommel Saddle Seat/1600 Simple Seat Back Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Tractor Seat with bracket <input type="checkbox"/> Large Pommel Saddle Seat <input type="checkbox"/> Saddle Seat		
<input type="checkbox"/> Option #3 – Alternate Seating System: Gray Bucket Seat		
1416 Accessories: <input type="checkbox"/> Heavy-Duty Push Bar <input type="checkbox"/> Pedal Toe Pulley <input type="checkbox"/> H-Harness - 11.5" <input type="checkbox"/> Rear Steering Kit <input type="checkbox"/> Calf & Leg Supports <input type="checkbox"/> ½" Exercise Pedals <input type="checkbox"/> Hi-Rise Handlebars: <input type="checkbox"/> 8" <input type="checkbox"/> 12" <input type="checkbox"/> Solid Tires <input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5" <input type="checkbox"/> 9" extender tube		

<input type="checkbox"/> 1420 ProSeries FOOT TRYKE (20" wheels, fixed or freewheel drive) – 50-FC-1420 Arm Length 20-30", Leg Length 29-35" NEW: Standard with Large Pommel Saddle Seat and safety flag		
<input type="checkbox"/> Option #1 – Standard Seating System: 1400 ProSeries Seat Back System (push grip, backpad, large pommel & 2 laterals) Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Tractor Seat w/ bracket <input type="checkbox"/> Medium Pommel Saddle Seat <input type="checkbox"/> Saddle Seat Trunk Support Options: <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Full Padded Back <input type="checkbox"/> ProSeries Lumbar Pad <input type="checkbox"/> Recumbent 10 Degree Seat Post		
<input type="checkbox"/> Option #2 – Alternate Seating System: Large Pommel Saddle Seat/1600 Simple Seatback Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Tractor Seat w/ bracket <input type="checkbox"/> Medium Pommel Saddle Seat <input type="checkbox"/> Saddle Seat		
<input type="checkbox"/> Option #3 – Alternate Seating System: Gray Bucket Seat		
1420 Accessories: <input type="checkbox"/> Heavy-Duty Push Bar <input type="checkbox"/> Pedal Toe Pulley <input type="checkbox"/> H-Harness - 11.5" <input type="checkbox"/> Rear Steering Kit <input type="checkbox"/> Calf & Leg Supports <input type="checkbox"/> ½" Exercise Pedals <input type="checkbox"/> 12" Hi-Rise Handlebars <input type="checkbox"/> Solid Tires <input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5" <input type="checkbox"/> 9" extender tube		

<input type="checkbox"/> 1420XL ProSeries FOOT TRYKE (20" wheels, fixed/freewheel drive) – 50-FC-1420-XL Arm Lgth 22-34", Leg Lgth 30-45" NEW: Standard with Large Pommel Saddle Seat and safety flag		
<input type="checkbox"/> Option #1 – Standard Seating System: 1400 ProSeries Seat Back System (push grip, backpad, large pommel & 2 laterals) Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Tractor Seat w/ bracket <input type="checkbox"/> Medium Pommel Saddle Seat <input type="checkbox"/> Saddle Seat Trunk Support Options: <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Full Padded Back <input type="checkbox"/> ProSeries Lumbar Pad <input type="checkbox"/> Recumbent 10 Degree Seat Post		
<input type="checkbox"/> Option #2 – Alternate Seating System: Large Pommel Saddle Seat/1600 Simple Seat Back Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Tractor Seat w/ bracket <input type="checkbox"/> Medium Pommel Saddle Seat <input type="checkbox"/> Saddle Seat		
<input type="checkbox"/> Option #3 – Alternate Seating System: Large Wheelchair Seat Wheelchair Seat Accessories: <input type="checkbox"/> Wheelchair Seat Bar Ends <input type="checkbox"/> Swing Away Arms		
1420 XL Accessories: <input type="checkbox"/> Heavy-Duty Push Bar <input type="checkbox"/> Toe Pedal Pulley <input type="checkbox"/> H-Harness - 11.5" <input type="checkbox"/> Rear Steering Kit <input type="checkbox"/> Calf & Leg Supports <input type="checkbox"/> ½" Exercise Pedals <input type="checkbox"/> Hi-Rise Handlebars <input type="checkbox"/> 12" <input type="checkbox"/> 23" <input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5" <input type="checkbox"/> 9" extender tube		

Before an order can be placed, we must receive a fully completed Request & Waiver Form, Assessment Form and Model Selection Form.

Please e-mail completed forms to caprocktrykes@gmail.com or mail copy to P.O. Box 93871, Lubbock, Tx 79423.

You will receive either a confirmation e-mail from caprocktrykes@gmail.com or a phone call from a member to inform you that your application has been received.

Recipient's Name: _____

FOOT TRYKES

<input type="checkbox"/> JT-2000 FOOT TRYKE (Recumbent 14-Speed) – 50-FC-2000 (advanced riders only) Arm Length 20-28", Leg Length 30-41"				
JT-2000 Accessories:	<input type="checkbox"/> Exercise Pedals	<input type="checkbox"/> XL Exercise Pedals	<input type="checkbox"/> Expanding Pedals	<input type="checkbox"/> Digital Speedometer
	<input type="checkbox"/> Dual Hand Brake	<input type="checkbox"/> Toe Clips	<input type="checkbox"/> Basket	<input type="checkbox"/> Rearview Mirror

<input type="checkbox"/> JT-2300-USS FOOT TRYKE (Recumbent 14-Speed with Under Seat Steering) – 50-FC-2300-USS (advanced riders only) Arm Length 20-28", Leg Length 30-41"				
JT-2300-USS Accessories:	<input type="checkbox"/> Exercise Pedals	<input type="checkbox"/> XL Exercise Pedals	<input type="checkbox"/> Expanding Pedals	<input type="checkbox"/> Digital Speedometer
	<input type="checkbox"/> Dual Hand Brake	<input type="checkbox"/> Toe Clips	<input type="checkbox"/> Basket	<input type="checkbox"/> Rearview Mirror

<input type="checkbox"/> TP-3000 TADPOLE FOOT TRYKE (Recumbent 21-Speed) – 50-FC-3000 (advanced riders only) Arm Length 17-29", Leg Length 36-45"				
Tadpole Accessories:	<input type="checkbox"/> Exercise Pedals	<input type="checkbox"/> XL Exercise Pedals	<input type="checkbox"/> Expanding Pedals	
	<input type="checkbox"/> Dual Hand Brake	<input type="checkbox"/> Toe Clips	<input type="checkbox"/> Rearview Mirror	

Therapist Assembly Notes or Comments (if any):

Before an order can be placed, we must receive a fully completed Request & Waiver Form, Assessment Form and Model Selection Form.
Please e-mail completed forms to caprocktrykes@gmail.com or mail copy to P.O. Box 93871, Lubbock, Tx 79423.
You will receive either a confirmation e-mail from caprocktrykes@gmail.com or a phone call from a member to inform you that your application has been received.

Recipient's Name:	
-------------------	--

HAND TRYKES

<input type="checkbox"/> AM-10 HAND TRYKE (12" wheels, 2.5" crank arms & permanent foot platform) 50-HFC-0105 Arm length 13-17" New: Now standard with Rear Steering (converts to push bar by inserting pin) and safety flag.							
<input type="checkbox"/> Option #1 – Standard Seating System: Blue Bucket Seat							
<input type="checkbox"/> Option #2 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: <table style="width: 100%; border: none;"> <tr> <td><input type="checkbox"/> Notched Seat Bottom Cushion*</td> <td><input type="checkbox"/> Laterals</td> <td><input type="checkbox"/> Head Rest</td> </tr> <tr> <td><input type="checkbox"/> Notched Wide Bottom Cushion*</td> <td><input type="checkbox"/> XL Seat Back**</td> <td></td> </tr> </table>		<input type="checkbox"/> Notched Seat Bottom Cushion*	<input type="checkbox"/> Laterals	<input type="checkbox"/> Head Rest	<input type="checkbox"/> Notched Wide Bottom Cushion*	<input type="checkbox"/> XL Seat Back**	
<input type="checkbox"/> Notched Seat Bottom Cushion*	<input type="checkbox"/> Laterals	<input type="checkbox"/> Head Rest					
<input type="checkbox"/> Notched Wide Bottom Cushion*	<input type="checkbox"/> XL Seat Back**						
AM-10 Accessories:	<input type="checkbox"/> Vertical Hand Grips <input type="checkbox"/> H-Harness - 11.5" <input type="checkbox"/> Heavy Duty Push Bar (<i>notched seats</i>) <input type="checkbox"/> Separator Cube						

<input type="checkbox"/> AM-12Small HAND TRYKE (12" wheels, 3" crank arms & permanent foot platform) 50-HFC-0110 Arm length 15-20" New: Now standard with Saddle Seat, Rear Steering (converts to push bar by inserting pin) and safety flag.							
<input type="checkbox"/> Option #1 – Standard Seating System: Saddle Seat & 1600 Simple Seat Back							
<input type="checkbox"/> Option #2 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: <table style="width: 100%; border: none;"> <tr> <td><input type="checkbox"/> Notched Seat Bottom Cushion*</td> <td><input type="checkbox"/> Laterals</td> <td><input type="checkbox"/> Head Rest</td> </tr> <tr> <td><input type="checkbox"/> Notched Wide Bottom Cushion*</td> <td><input type="checkbox"/> XL Seat Back**</td> <td></td> </tr> </table>		<input type="checkbox"/> Notched Seat Bottom Cushion*	<input type="checkbox"/> Laterals	<input type="checkbox"/> Head Rest	<input type="checkbox"/> Notched Wide Bottom Cushion*	<input type="checkbox"/> XL Seat Back**	
<input type="checkbox"/> Notched Seat Bottom Cushion*	<input type="checkbox"/> Laterals	<input type="checkbox"/> Head Rest					
<input type="checkbox"/> Notched Wide Bottom Cushion*	<input type="checkbox"/> XL Seat Back**						
<input type="checkbox"/> Option #3 – Alternate Seating System: Blue Bucket Seat							
AM-12S Accessories:	<input type="checkbox"/> Vertical Hand Grips <input type="checkbox"/> H-Harness - 11.5" <input type="checkbox"/> Heavy Duty Push Bar (<i>notched seats</i>) <input type="checkbox"/> Separator Cube						

<input type="checkbox"/> AM-12 HAND TRYKE – (12" wheels, 4" crank arms & permanent foot platform) 50-HC-0210 Arm Length 14-23" New: Now standard with Rear Steering (converts to push bar by inserting pin) and safety flag.								
<input type="checkbox"/> Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seat Back Seat Bottom Alternates: <table style="width: 100%; border: none;"> <tr> <td><input type="checkbox"/> Bench Seat</td> <td><input type="checkbox"/> Pommel Saddle Seat</td> <td><input type="checkbox"/> Medium (11.5x10")</td> <td><input type="checkbox"/> Large(14x13")</td> </tr> </table>		<input type="checkbox"/> Bench Seat	<input type="checkbox"/> Pommel Saddle Seat	<input type="checkbox"/> Medium (11.5x10")	<input type="checkbox"/> Large(14x13")			
<input type="checkbox"/> Bench Seat	<input type="checkbox"/> Pommel Saddle Seat	<input type="checkbox"/> Medium (11.5x10")	<input type="checkbox"/> Large(14x13")					
<input type="checkbox"/> Option #2 – Alternate Seating System: 1400 ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals) Seat Bottom Alternates: <table style="width: 100%; border: none;"> <tr> <td><input type="checkbox"/> Bench Seat</td> <td><input type="checkbox"/> Pommel Saddle Seat</td> <td><input type="checkbox"/> Medium (11.5x10")</td> <td><input type="checkbox"/> Large(14x13")</td> </tr> </table> Trunk Support Options: <table style="width: 100%; border: none;"> <tr> <td><input type="checkbox"/> ProSeries Full Padded Back</td> <td><input type="checkbox"/> ProSeries Head Rest</td> <td><input type="checkbox"/> ProSeries Lumbar Pad</td> </tr> </table>		<input type="checkbox"/> Bench Seat	<input type="checkbox"/> Pommel Saddle Seat	<input type="checkbox"/> Medium (11.5x10")	<input type="checkbox"/> Large(14x13")	<input type="checkbox"/> ProSeries Full Padded Back	<input type="checkbox"/> ProSeries Head Rest	<input type="checkbox"/> ProSeries Lumbar Pad
<input type="checkbox"/> Bench Seat	<input type="checkbox"/> Pommel Saddle Seat	<input type="checkbox"/> Medium (11.5x10")	<input type="checkbox"/> Large(14x13")					
<input type="checkbox"/> ProSeries Full Padded Back	<input type="checkbox"/> ProSeries Head Rest	<input type="checkbox"/> ProSeries Lumbar Pad						
<input type="checkbox"/> Option #3 – Alternate Seating System: Gray Bucket Seat								
<input type="checkbox"/> Option #4 – Alternate Seating System: Snappy Seat System Snappy Seat Accessories: <table style="width: 100%; border: none;"> <tr> <td><input type="checkbox"/> Notched Seat Bottom Cushion*</td> <td><input type="checkbox"/> Laterals</td> <td><input type="checkbox"/> Head Rest</td> </tr> <tr> <td><input type="checkbox"/> Notched Wide Bottom Cushion*</td> <td><input type="checkbox"/> XL Seat Back**</td> <td><input type="checkbox"/> Separator Cube</td> </tr> </table>		<input type="checkbox"/> Notched Seat Bottom Cushion*	<input type="checkbox"/> Laterals	<input type="checkbox"/> Head Rest	<input type="checkbox"/> Notched Wide Bottom Cushion*	<input type="checkbox"/> XL Seat Back**	<input type="checkbox"/> Separator Cube	
<input type="checkbox"/> Notched Seat Bottom Cushion*	<input type="checkbox"/> Laterals	<input type="checkbox"/> Head Rest						
<input type="checkbox"/> Notched Wide Bottom Cushion*	<input type="checkbox"/> XL Seat Back**	<input type="checkbox"/> Separator Cube						
AM-12 Accessories:	<input type="checkbox"/> Vertical Hand Grips <input type="checkbox"/> H-Harness - 11.5" <input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5" <input type="checkbox"/> 9" Extender Tube							

*cannot be used in conjunction with rear steering. **cannot be used in conjunction with an h-harness, head rest or laterals.

Before an order can be placed, we must receive a fully completed Request & Waiver Form, Assessment Form and Model Selection Form.

Please e-mail completed forms to caprocktrykes@gmail.com or mail copy to P.O. Box 93871, Lubbock, Tx 79423.

You will receive either a confirmation e-mail from caprocktrykes@gmail.com or a phone call from a member to inform you that your application has been received.

<input type="checkbox"/> AM-16 HAND TRYKE (16" wheels, 5" crank arms & permanent foot platform) 50-HFC-0411 Arm length 18-27" New: Now standard with Rear Steering (converts to push bar by inserting pin) and safety flag.	
<input type="checkbox"/> Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seat Back Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13")	
<input type="checkbox"/> Option #2 – Upgraded Seating System: 1400 ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals) Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Padded Tractor Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13") Trunk Support Options: <input type="checkbox"/> ProSeries Full Padded Back <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Lumbar Pad	
<input type="checkbox"/> Option #3 – Upgraded Seating System: Gray Bucket Seat	
AM-16 Accessories: <input type="checkbox"/> Vertical Hand Grips <input type="checkbox"/> H-Harness - 11.5" <input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5" <input type="checkbox"/> 9" Extender Tube	

Recipient's Name:

HAND TRYKES

<input type="checkbox"/> NEW: 1020 "Junior" HAND TRYKE (20" wheels & 3-Speed drive train) 50-HC-1020 Arm length 19-30"	
<input type="checkbox"/> Option #1 – Standard Seating: Small Wheelchair Seat (3.5" narrower in width and depth than the Large) Seat Alternates: <input type="checkbox"/> Large Wheelchair Seat	
1020 Accessories: <input type="checkbox"/> Wheelchair Seat Bar Ends <input type="checkbox"/> Swing-Away Arms <input type="checkbox"/> Rearview Mirror	

<input type="checkbox"/> 1024 HAND TRYKE – (24" wheels & 3-Speed drive train) 50-HC-1024 Arm length 22-26"	
<input type="checkbox"/> Option #1 – Standard Seating: Large Wheelchair Seat (3.5" broader in width and depth than the Small) Seat Alternates: <input type="checkbox"/> Small Wheelchair Seat	
1024 Accessories: <input type="checkbox"/> Wheelchair Seat Bar Ends <input type="checkbox"/> Swing-Away Arms <input type="checkbox"/> Rearview Mirror	

Therapist Assembly Notes or Comments (if any):
--

Before an order can be placed, we must receive a fully completed Request & Waiver Form, Assessment Form and Model Selection Form.

Please e-mail completed forms to caprocktrykes@gmail.com or mail copy to P.O. Box 93871, Lubbock, Tx 79423.

You will receive either a confirmation e-mail from caprocktrykes@gmail.com or a phone call from a member to inform you that your application has been received.

Recipient's Name: _____

HAND-FOOT TRYKES

<input type="checkbox"/> AM-10 HAND-FOOT TRYKE (10" wheels & 2.5" crank arms) 50-HFC-0105 Arm length 13-17", leg length 15-21" New: Now standard with Rear Steering (converts to push bar by inserting pin) and safety flag.			
<input type="checkbox"/> Option #1 – Standard Seating System: Blue Bucket Seat			
<input type="checkbox"/> Option #2 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest <input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back**			
AM-10 Accessories:	<input type="checkbox"/> AM Toe Pulley	<input type="checkbox"/> H-Harness - 11.5"	
	<input type="checkbox"/> Heavy-Duty Push Bar (<i>notched seats</i>)	<input type="checkbox"/> Vertical Hand Grips	<input type="checkbox"/> Separator Cube

<input type="checkbox"/> AM-12Small HAND-FOOT TRYKE (12" wheels & 3" crank arms) 50-HFC-0110 Arm length 15-20", leg length 19-24" New: Now standard with Saddle Seat, Rear Steering (converts to push bar by inserting pin) and safety flag.			
<input type="checkbox"/> Option #1 – Standard Seating System: Saddle Seat & 1600 Simple Seat Back			
<input type="checkbox"/> Option #2 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest <input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back**			
<input type="checkbox"/> Option #3 – Upgraded Seating System: Blue Bucket Seat			
AM-12S Accessories:	<input type="checkbox"/> AM Pedal Toe Pulley	<input type="checkbox"/> H-Harness - 11.5"	
	<input type="checkbox"/> Heavy-Duty Push Bar (<i>notched seats</i>)	<input type="checkbox"/> Vertical Hand Grips	<input type="checkbox"/> Separator Cube
		<input type="checkbox"/> 9/16" Exercise pedals	<input type="checkbox"/> Solid Rear Tires

*cannot be used in conjunction with rear steering. **cannot be used in conjunction with an h-harness, head rest or laterals.

<input type="checkbox"/> AM-12 HAND-FOOT TRYKE (12" wheels & 4" crank arms) 50-HFC-0210 Arm length 14-23", leg length 21-29" New: Now standard with Rear Steering (converts to push bar by inserting pin) and safety flag.			
<input type="checkbox"/> Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seat Back Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13")			
<input type="checkbox"/> Option #2 – Upgraded Seating System: ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals) Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13") Trunk Support Upgrades: <input type="checkbox"/> ProSeries Full Padded Back <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Lumbar Pad			
<input type="checkbox"/> Option #3 – Upgraded Seating System: Gray Bucket Seat			
<input type="checkbox"/> Option #4 – Upgraded Seating System: Snappy Seat System (<i>too small for many AM-12 riders; check seat measurements.</i>) Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest <input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back** <input type="checkbox"/> Separator Cube			
AM-12 Accessories:	<input type="checkbox"/> AM Pedal Toe Pulley	<input type="checkbox"/> H-Harness - 11.5"	
	<input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5"	<input type="checkbox"/> 9" Extender Tube	<input type="checkbox"/> Vertical Hand Grips
		<input type="checkbox"/> 9/16" Exercise pedals	<input type="checkbox"/> Solid Rear Tires

Before an order can be placed, we must receive a fully completed Request & Waiver Form, Assessment Form and Model Selection Form.

Please e-mail completed forms to caprocktrykes@gmail.com or mail copy to P.O. Box 93871, Lubbock, Tx 79423.

You will receive either a confirmation e-mail from caprocktrykes@gmail.com or a phone call from a member to inform you that your application has been received.

Recipient's Name:

AM-16 HAND-FOOT TRYKE (16" wheels & 5" crank arms) **50-HFC-0411** Arm length 18-27", leg length 24-36"
New: Now standard with Rear Steering (converts to push bar by inserting pin) and safety flag.

Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seat Back
Seat Bottom Alternates: Bench Seat Tractor Seat Pommel Saddle Seat Medium (11.5x10") Large(14x13")

Option #2 – Upgraded Seating System: 1400 ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals)
Seat Bottom Upgrades: Bench Seat Tractor Seat Pommel Saddle Seat Medium (11.5x10") Large(14x13")
Trunk Support Upgrades: ProSeries Full Padded Back ProSeries Head Rest ProSeries Lumbar Pad

Option #3 – Upgraded Seating System: Gray Bucket Seat

AM-16 Accessories: AM Pedal Toe Pulley H-Harness - 11.5"
 Knee Separator: 2" 5" 9" Extender Tube XL Exercise pedals
 Solid Rear Tires 9/16" Exercise pedals Vertical Hand Grips

AM-20 HAND-FOOT TRYKE (Recumbent Single Speed) – **50-HFC-0610** Arm length 18-30", leg length 32-46"
Tryke is currently being used ONLY for research.

AM-20 Accessories: Exercise Pedals XL Exercise Pedals Expanding Pedals
 Wheelchair Seat Bar Ends Swing Away Arms Calf & Leg Supports (Expanding pedals)

Therapist Assembly Notes or Comments (if any):

Before an order can be placed, we must receive a fully completed Request & Waiver Form, Assessment Form and Model Selection Form.
Please e-mail completed forms to caprocktrykes@gmail.com or mail copy to P.O. Box 93871, Lubbock, Tx 79423.
You will receive either a confirmation e-mail from caprocktrykes@gmail.com or a phone call from a member to inform you that your application has been received.